


Ultra Pure Human Apolipoprotein AI (Apo AI)

Product Specification

Catalog No:	11P-UP210, 1 mg
Protein Concentration:	1mg / ml, determined by the Lowry method
Source:	From fresh human plasma that has tested negative for Hepatitis C, HIV-I and HIV-II antibodies as well as Hepatitis surface antigens.
Purification:	After series ultracentrifugations, High Density Lipoprotein (HDL) is isolated from human plasma. Apo AI is purified from delipidated HDL, followed by gel-filtration.
Purity:	> 99 % by SDS-PAGE
Preservation:	20 mM Tris-HCl, 0.14 M NaCl, 0.02 % NaN ₃ , 0.5 mM EDTA, pH 7.4.
Storage:	-20° C for long-term storage, 2-8° C for short-term storage. Aliquot to avoid repeated freezing and thawing.



**The products are for research or manufacturing use only, not for use in human therapeutic or diagnostic applications.*

