



## Nitrated Bovine Serum Albumin (Nitrated BSA), with 50% glycerol

# Product Specification

Catalog No:	20P-NY-BS105, 0.5 mg
Protein Concentration:	1 mg / ml, determined by the Lowry method
Source:	Bovine Serum Albumin (purchased from Boehringer Mannheim) modified with nitrotyrosine.
Preservation:	10 mM PBS, 0.15 M NaCl 0.5 mM EDTA, 0.02 % NaN <sub>3</sub> , pH 7.4. In addition, preserved with 50 % glycerol.
Storage:	-20° C for short and long-term storage. Aliquot to avoid repeated freezing and thawing.

*\*The products are for research or manufacturing use only, not for use in human therapeutic or diagnostic applications.*

